

# A New Enzymatic Cycling Method for Determination of Total Bile Acids in Serum

H. Mueller<sup>1</sup>, I. Delseith<sup>1</sup>, A. Klich<sup>1</sup>, T. Leber<sup>1</sup>, S. Runkel<sup>1</sup>, M. Grimmmler<sup>1,2</sup>

<sup>1</sup>DiaSys Diagnostic Systems GmbH, Alte Strasse 9, 65558 Holzheim, Germany

<sup>2</sup>Fresenius University of Applied Science, Limburger Strasse 2, 65510 Idstein, Germany

## BACKGROUND

Serum total bile acids are a sensitive indicator of liver function and can be used for the diagnosis and prognosis of various liver diseases such as viral hepatitis, intrahepatic cholestasis of pregnancy, liver cancer or toxic hepatic damage. Commercially available assays show limitations regarding the detection of various bile acids. The objective of this study was to compare the new DiaSys total bile acids test with other assays in the market, in particular with regard to the detection of diagnostically relevant primary and secondary bile acids and their conjugates.

## METHODOLOGY

The kinetic DiaSys total bile acids test, Bile Acids 21 FS, is based on a cycling reaction which reversibly changes bile acids to an oxidized form, catalyzed by a specific 3- $\alpha$ -hydroxysteroid-dehydrogenase. In this reaction Thio-NAD is reduced to Thio-NADH. In a second reaction step, oxidized bile acids are reduced by the same enzyme with subsequent reduction of NADH to NAD. The rate of Thio-NADH formation is determined by measuring a specific change of absorbance at 405/415 nm, which is directly proportional to the concentration of bile acids in the sample material. Comparative studies were performed on Hitachi 917 with 52 serum samples (mix of native and spiked samples; spiking was done with equal amounts of aqueous glycocholic acid, glycochenodeoxycholic acid, taurocholic acid and taurodesoxycholic acid solutions) [Fig. 1,2]. Data have been evaluated by using regression analysis according to Passing and Bablok. Recovery studies were performed using 50  $\mu$ M aqueous bile acids solutions [Table 1]. In-series precision studies were performed according to CLSI protocol (EP5-A3) [Table 2].

## RESULTS

Bile Acids 21 FS is a liquid stable, 2-component reagent with a sample reagent ratio of 4:270:90  $\mu$ L (on Hitachi 917). Linearity is given up to 220  $\mu$ mol/L. The reagent shows good correlation for sera in comparison studies with commercially available methods [Fig. 1,2]. The DiaSys test outperformed the competitors' assays with respect to the recovery (Mean of deviations to 100% = 5.10% vs. 8.12% (competitor A) vs. 23.67% (competitor B)) [Table 1]. Furthermore, Bile Acids 21 FS shows a good in-series precision with a CV of  $\leq$  2.89% [Table 2].

| Recovery in aqueous bile acids solution [50 $\mu$ M] | DiaSys Cycling Method |       | Competitor A Cycling Method |       | Competitor B Cycling Method |       |
|--|-----------------------|-------|-----------------------------|-------|-----------------------------|-------|
|  | [ $\mu$ mol/L]        | [%]   | [ $\mu$ mol/L]              | [%]   | [ $\mu$ mol/L]              | [%]   |
| Glycocholic acid                                     | 50.3                  | 100.5 | 49.0                        | 97.9  | 33.8                        | 67.5  |
| Glycochenodeoxycholic acid                           | 46.0                  | 92.0  | 52.0                        | 104.0 | 50.9                        | 101.9 |
| Taurochenodesoxycholic acid                          | 48.0                  | 96.1  | 42.1                        | 84.3  | 44.8                        | 89.6  |
| Taurocholic acid                                     | 49.6                  | 99.2  | 44.8                        | 89.7  | 31.9                        | 63.7  |
| Chenodeoxycholic acid                                | 52.0                  | 104.0 | 55.8                        | 111.7 | 54.9                        | 109.8 |
| Cholic acid  | 53.6                  | 107.3 | 50.4                        | 100.9 | 35.8                        | 71.5  |
| Desoxycholic acid                                    | 53.9                  | 107.8 | 55.8                        | 111.7 | 69.7                        | 139.4 |
| Taurodesoxycholic acid                               | 51.4                  | 102.7 | 48.0                        | 96.0  | 59.7                        | 119.5 |
| Lithocholic acid                                     | 44.6                  | 89.1  | 43.7                        | 87.3  | 32.7                        | 65.4  |

Table 1: Recovery of primary and secondary bile acids




Figure 1: Competitor A - comparison total bile acids




Figure 2: Competitor B - comparison total bile acids

| Mean [ $\mu$ mol/L] | SD [ $\mu$ mol/L] | CV [%] |
|---------------------|-------------------|--------|
| 7.37                | 0.21              | 2.89   |
| 16.2                | 0.27              | 1.67   |
| 36.7                | 0.42              | 1.13   |

Table 2: In-series precision (n=20)

## CONCLUSIONS

All results clearly show that Bile Acids 21 FS meets the requirements of a sensitive test to analyze all stages of impaired liver function. Diagnostically relevant primary and secondary bile acids as well as their conjugates are reliably recovered. The performance of the test is highly competitive in comparison to competitor products.

## REFERENCES

- Zhang GH *et al.* (2005). An enzymatic cycling method for the determination of serum total bile acids with recombinant 3 $\alpha$ -hydroxysteroid dehydrogenase. *Biochem Biophys Res Commun*; 326: 87–92.
- Hofmann AF (1999). The continuing importance of bile acids in liver and intestinal disease. *Arch Intern Med*; 159:2647–2658.
- Thomas L ed. *Clinical Laboratory Diagnostics*. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 1998:1341-44.
- Neale G *et al.* (1971). Serum bile acids in liver disease. *Gut*; 12: 145-152.