

#### Leading Technology in Fluid-Stable Reagents from DiaSys


- More than 25 years experience in development and production of clinical chemistry tests
- Premium service in technics, applications and after sales
- Quality products made in Germany
- High performance, ready-to-use reagents with minimized interferences, long shelf life and onboard stability as well as traceability to international references
- Perfectly matched fluid-stable reagents, calibrators and controls
- High grade raw materials from traceable origin
- Processes and resources certified according to ISO 13485, fulfilling highest quality standards
- Sustainable processes and products preserve the environment

DiaSys offers reagent kits for manual and automated use plus appropriate calibrators and controls. Detailed information about  $\alpha$ -amylase, pancreatic amylase and lipase test is available on our website [www.diasys-diagnostics.com/products/reagents](http://www.diasys-diagnostics.com/products/reagents) and in our product catalog.

Handed over by:



Climate neutral print  
(Carbon neutral)  
ID: DE-626-983465  
Printed on FSC certified paper



**DiaSys**  
**Diagnostic Systems GmbH**  
Alte Strasse 9  
65558 Holzheim  
Germany

Phone: +49 6432 9146-0  
Fax: +49 6432 9146-32  
E-Mail: [reagents@diasys.de](mailto:reagents@diasys.de)  
[www.diasys-diagnostics.com](http://www.diasys-diagnostics.com)




820020 | September 2019

CHOOSING QUALITY.

# Pancreatic Diagnostics

## Fast. Easy. Reliable.



- Lipase
- $\alpha$ -Amylase
- Pancreatic Amylase



CHOOSING QUALITY.

# Clinical Use

Pancreatic disorders, both inflammatory and neoplastic, are increasingly occurring in industrialized countries. Clinical chemistry tests for α-amylase, pancreatic amylase and lipase are used in detection and follow-up of acute pancreatitis and diagnosis of abdominal disorders. Amylase in serum consists of two isoenzymes: Salivary amylase originating from the salivary glands and pancreatic amylase produced by the pancreas, of which only pancreatic amylase is specific for pancreatic disorders. Lipase represents a further specific parameter for pancreatic diseases. Due to its long half-life, lipase values are elevated longer after a pancreatitis than those of alpha and pancreatic amylase. The combination of both, lipase and pancreatic amylase, highly increases the sensitivity for the detection of pancreatitis.

# Methods

**Lipase DC FS:** The DiaSys test is a state-of-the-art direct colorimetric assay. It uses a substrate in a microemulsion with colipase and bile acids, thus simulating physiological conditions and, therefore, avoiding non-specific enzyme reactions.

**α-Amylase CC FS and Pancreatic amylase CC FS:** Both tests have the same reaction principle. The DiaSys method uses an Ethylidene protected p-Nitrophenyl-α-D-maltoheptaoside (EPS-G7) as substrate and a bacterial α-glucosidase as recommended by the IFCC (International Federation of Clinical Chemistry). This combination leads to a very stable reagent with fast and complete color release. For the specific measurement of pancreatic amylase, two monoclonal antibodies are used to inhibit the salivary amylase.

**DiaSys Reagents for Diagnostics of Pancreatic Diseases**

- Liquid stable up to the expiry date, even after opening
- Ready to use – no reconstitution, directly available for stat analysis
- Adaptable to automated clinical chemistry analyzers, validated applications available

**Lipase DC FS**

- Direct colorimetric method: High signal, high measuring range
- Specific substrate in microemulsion with bile acids and colipase avoids unspecific reactions
- Use of serum and heparin plasma
- Linearity up to 300 U/L\*
- Calibration/onboard stability 12 weeks (BioMajesty® JCA-BM6010/C)

**α-Amylase CC FS**

- IFCC recommended EPS-G7 method
- Linearity up to 2000 U/L\*
- Good correlation to other methods (e.g. CNP-G3)
- Use of serum, plasma and urine
- Mixed reagent stability 6 months
- Calibration/onboard stability 6 weeks (BioMajesty® JCA-BM6010/C)


**Pancreatic amylase CC FS**

- EPS-G7 method, based on IFCC recommendation
- Immunoinhibition of salivary amylase: Specific for pancreatic diagnostics
- Use of serum, plasma and urine
- Linearity up to 2000 U/L\*
- Calibration/onboard stability 6 weeks (BioMajesty® JCA-BM6010/C)

\* May vary according to specific applications.


**Lipase DC FS:** Precision studies and method comparison to competitor Lipase (on BioMajesty® JCA-BM6010/C)

Precision			
Intra-assay	Mean (U/L)	Mean (μkat/L)	CV (%)
Sample 1	30.9	0.52	1.26
Sample 2	60.9	1.02	0.611
Sample 3	286	4.78	0.263
Total	Mean (U/L)	Mean (μkat/L)	CV (%)
Sample 1	30.2	0.50	2.01
Sample 2	59.9	1.00	1.20
Sample 3	284	4.74	1.10



**α-Amylase CC FS:** Precision studies and method comparison to competitor α-Amylase (on BioMajesty® JCA-BM6010/C)

Precision			
Intra-assay	Mean (U/L)	Mean (μkat/L)	CV (%)
Sample 1	36.9	0.62	1.86
Sample 2	74.9	1.25	1.11
Sample 3	238	3.97	0.78
Inter-assay	Mean (U/L)	Mean (μkat/L)	CV (%)
Sample 1	37.3	0.62	1.83
Sample 2	119	1.98	1.59
Sample 3	242	4.04	1.90



**Pancreatic amylase CC FS:** Precision studies and instrument comparison of DiaSys Pancreatic amylase CC FS on Hitachi 917 vs. BioMajesty® JCA-BM6010/C

Precision			
Intra-assay	Mean (U/L)	Mean (μkat/L)	CV (%)
Sample 1	26.3	0.44	0.90
Sample 2	115	1.91	0.64
Sample 3	252	4.20	0.61
Inter-assay	Mean (U/L)	Mean (μkat/L)	CV (%)
Sample 1	44.9	0.75	1.40
Sample 2	142	2.37	0.86
Sample 3	247	4.13	0.59

