

# Bicarbonate FS

## Keep the Balance



- Fluid-Stable, Ready-to-Use Monoreagent
- Superior On-Board and Calibration Stability
- Excellent Performance

**DiaSys**

DiaSys  
Diagnostic Systems

CHOOSING QUALITY.

# Clinical Relevance

Bicarbonate, also termed hydrogencarbonate, is a salt derived from carbonic acid. Since the anions of this salt ( $\text{HCO}_3^-$  /  $\text{CO}_3^{2-}$ ) can act as a proton donor or acceptor, they fulfill the essential physiological function as a buffer

system to ensure the maintenance of the acid-base balance in blood. Deviating values indicate disorders associated with disturbances of metabolic and respiratory systems like acidosis or alkalosis.

## Reaction Principle

Phosphoenolpyruvate +  $\text{HCO}_3^-$

Phosphoenolpyruvate Carboxylase +  $\text{Mg}^{2+}$

Oxaloacetate +  $\text{H}_2\text{PO}_4^-$

Oxaloacetate + NADH analog red.

Malate Dehydrogenase

Malate + NADH analog

The decrease in absorbance measured at 405/415 nm is directly proportional to the bicarbonate concentration in the sample.

Intra-assay precision n = 20	Mean (mmol/L)	SD (mmol/L)	CV (%)
Sample 1	17.6	0.14	0.80
Sample 2	19.9	0.16	0.80
Sample 3	30.1	0.28	0.93

Inter-assay precision n = 20	Mean (mmol/L)	SD (mmol/L)	CV (%)
Sample 1	16.8	0.53	3.16
Sample 2	20.3	0.49	2.40
Sample 3	30.0	0.68	2.26

## Excellent Performance


- Synthetic NADH analog ensures excellent stability
- 3 weeks on-board and calibration stability
- No interferences by ascorbate up to 30 mg/dL, bilirubin up to 50 mg/dL, hemoglobin up to 500 mg/dL and lipemia up to 1400 mg/dL triglycerides
- Wide measuring range 4 – 50 mmol/L
- Easy applicable to various clinical chemistry analyzers
- Fluid-stable, ready-to-use monoreagent
- Very good comparison versus competitor test



Climate-neutral print  
(Carbon neutral)

ID: DE 626 955 906

Printed on FSC-certified paper



**DiaSys**  
**Diagnostic Systems GmbH**  
Alte Strasse 9  
65558 Holzheim  
Germany

Phone: +49 (0) 64 32 / 91 46-0  
Fax: +49 (0) 64 32 / 91 46-32  
E-Mail: [reagents@diasys.de](mailto:reagents@diasys.de)  
[www.diasys-diagnostics.com](http://www.diasys-diagnostics.com)



CHOOSING QUALITY.