












## REAGENT KITS SUITABLE FOR VARIOUS CLINICAL CHEMISTRY ANALYZERS

Details on reagent bottles and barcodes as follows:

System Kits			
<p>Last 3 digits of Cat. No.</p> <hr/> <p>... 917</p>		<p>R1</p>  <p>R2</p> 	<p><b>Bottles fit for:</b></p> <p>Hitachi 917 Modular P Beckman Coulter AU series ADVIA 1650/1800/2400 Abbott ARCHITECT series Dirui CS series</p> <p><b>Barcodes fit for:</b></p> <p>Hitachi 917 Modular P</p>
<p>Last 3 digits of Cat. No.</p> <hr/> <p>...704</p>		<p>R1</p>  <p>R2</p> 	<p><b>Bottles fit for:</b></p> <p>Hitachi 704/717/902/911/912</p> <p><b>Barcodes fit for:</b></p> <p>Hitachi 911/912</p>
<p>Last 3 digits of Cat. No.</p> <hr/> <p>...930/931/ 935/936/ 937</p>		<p>R1 + R2</p> 	<p><b>Bottles fit for:</b></p> <p>Hitachi 704/717/902/911/912/917 Modular P Beckman Coulter AU series ADVIA 1650/1800/2400 Abbott ARCHITECT series Dirui CS series</p> <p><b>Barcodes fit for:</b></p> <p>Hitachi 911/912/917 Modular P</p>
<p>Last 3 digits of Cat. No.</p> <hr/> <p>...717</p>		<p>R1</p>  <p>R2</p> 	<p><b>Bottles fit for:</b></p> <p>Hitachi 717/911/912</p> <p><b>Barcodes fit for:</b></p> <p>Hitachi 911/912</p>